

Finest Swiss Made solutions... worldwide!



**POLYTRON®**  
DISPERSING AND MIXING TECHNOLOGY BY KINEMATICA

**PT 10-35 GT**

*THE new standard  
benchtop homogenizer  
with space saving slim line  
design, new easy-to-use  
stand and LED-Display*

*Safe and reliable work  
in the Lab*

Work vol.	ml	0.1 - 4000
Speed	min <sup>-1</sup>	0-30'000
Tip Speed	m/s	up to 40
Power	W	1200 / 800
Noise level	dB(A)	< 75
Display	-	LED
Regulation / Control	-	integrated digital speed controller with overheating, overload and blocking protection
Disp.Aggr.	Ø mm	5 - 36
Coupling	-	Type F Quick coupling, compatible with Type A and Type B

**UPGRADED**



**KINEMATICA AG**  
Dispensing and Mixing Technology



**KINEMATICA, INC.**  
Dispensing and Mixing Technology

[www.kinematica.ch](http://www.kinematica.ch)  
[www.kinematica-inc.com](http://www.kinematica-inc.com)